




MaloBacti Activation Guide

IMPORTANT!




MaloBacti strains are not build-up cultures. We deliver you the full cell count so all you have to do is activate the bacteria in water; that's why wine is not required for the activation. As always, when dealing with yeasts or bacteria it is extremely important to strictly follow the activation instructions.




STEP 1

PLEASE FOLLOW INSTRUCTIONS CAREFULLY.




STEP 2

Check wine parameters*




STEP 3

Add chamber **1** to 1L of non-chlorinated, non-distilled water at 23-28°C. Stir well. Then add chamber **2** stir well for 5-8 min.




STEP 4

Measure and record the pH of the suspension.




STEP 5

Activate the suspension in a warm environment (23-28°C) for 6-8hrs.



STEP 6

During activation stir suspension twice.



STEP 7

Within 6-8 hours the pH should drop to 3.6-3.8, the bacteria are now completely activated. Verify with a pH meter.



STEP 8

Stir suspension and inoculate directly into 2500L of wine. Mix well.



STEP 9

Maintain wine temperature between 18-20°C.

25hL pouch = 1L non-chlorinated water

250hL pouch = 10L non-chlorinated water

100kL pouch = 40L non-chlorinated water

After complete activation the suspension can be stored at 4-6°C for maximum of 5 days. Bring suspension to >17°C and stir before inoculation.

**Even if wine parameters are within range we can not guarantee MLF due to other inhibitory factors that can be present in all wines.*



Kauri Winemaking Innovation

Aus Tel: 1800 127 611

NZ Tel: 0800 528 749

Email: info@kauriwine.com

Web: www.kauriwine.com